VarseoWax Model

Seite / page

de	Gebrauchsanweisung	2
en	Instructions for use	3
fr	Notice d'utilisation	4
es	Instrucciones de uso	5
it	Istruzioni per l'uso	6
ru	Инструкция по применению	7

BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG
Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany
Tel. +49 421 20 28-0 · Fax +49 421 20 28-100
E-Mail: info@bego.com · www.bego.com

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VarseoWax Model

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Instructions for use

VarseoWax Model

This device is a resin for the 3D printing of dental models.

1. Device description

VarseoWax Model is an acrylic acid-based monomer for the manufacture of 3D printed dental models. The device can be used in printers with a wavelength range of 385 to 405 nm. The verified, compatible 3D print system components can be found on our website <https://www.bego.com/3d-printing/compatibility-overview/>

2. Restrictions of use

VarseoWax Model contains ingredients that could trigger allergic reactions. In cases of doubt, the allergy should be clarified and ruled out through a specific test prior to the application of this device.

3. Safety instructions

This device is manufactured and tested in accordance with the most stringent quality standards. In order to ensure optimum further processing, please read the information contained in these instructions for use carefully.

Improper use and failure to observe the information can have a detrimental effect on the quality.

The safety instructions and precautions contained in the instructions for use and safety data sheet for this device shall apply to the handling of liquid resin and printed objects that have not been post-cured (objects in "green condition").

For professional use only.

4. Side-effects and precautions

Precautions / Protection

Protective clothing must be worn when handling this device.

Safety goggles and nitrile gloves must be used. Further information on handling the device can be found in the safety data sheet and also downloaded from the BEGO Download Centre at www.bego.com. However, we cannot completely rule out the possibility that personal reactions (e.g. intolerances or allergic reactions) to individual components could occur in isolated cases. In such cases, the respective user should discontinue use of the device.

5. General information on handling

Delivery

VarseoWax Model is supplied in lightproof, sealed bottles.

Filling quantity

- REF 41140 = 1 kg
- REF 41141 = 250 g

Storage

This device must be stored in the original sealed bottle at room temperature (approx. 22°C) in a dark, dry place. It must be ensured that the temperature does not drop below +5°C and does not exceed +35°C!

Note: The minimum shelf life date printed on the device must be observed. Perfect processing cannot be guaranteed if devices which have exceeded their minimum shelf life date are used.

When not being used or when being stored, the printed models that have not been post-cured (objects in "green condition") must be stored away from light to prevent changes in form! The fully cured print objects must be stored at room temperature and protected from sources of light.

6. Processing requirements

a. Design

- Create the object (STL data set) with a commercial CAD software intended for dental applications.
- During the design process, observe the required minimum wall thicknesses for completed restorations.

Note: Ensure a minimum wall thickness of 2.5 mm for hollow print models.

b. Compatible 3D printers

In principle, VarseoWax Model can be processed on 3D print systems that use a wavelength of between 385 and 405 nm.

The verified, compatible 3D print system components can be found on our website <https://www.bego.com/3d-printing/compatibility-overview/>

c. Nesting and preparation for printing

- STL file import
- Manual/automatic rotation and placement
- Optimal alignment: Position the print object horizontally on the build platform
- Manual/automatic generation of support structures.

In order to ensure optimal processing, please read carefully the information contained in the instructions for use for the nesting software.

d. Tools, equipment and materials required for post-processing

- Stainless steel spatula
- Cleaning equipment (e.g. unheated ultrasonic bath)
- Cleaning solution (e.g. ethanol solution 96 % / isopropanol 99 %)
- Spray bottle for the cleaning solution
- Cutting disc or side cutters (to remove the support structure)

7. Processing

VarseoWax Model is compatible and can be used together with the 3D print system components of different suppliers. These compatible system components can be found on our website <https://www.bego.com/3d-printing/compatibility-overview/>

Note: All processing methods listed in these instructions for use serve as an example.

Printing

Follow the instructions in the operating manual for the 3D printer when preparing for the printing process. The ideal working temperature for VarseoWax Model is in the range between 18°C and 28°C.

Before use, the resin must have a homogeneous consistency. Before first use, shake the bottle well for approx. 2 min. When decanting, make sure that the printing resin is exposed to daylight for as short a period of time as possible. Mix the resin in the cartridge/resin tub if a transparent layer is visible on the surface. During the printing process, follow the operating instructions of the respective printer for further processing – selecting the resin, material parameters, setting up the print job.

Cleaning

Upon completion of printing, the print objects are released from the build platform using a spatula. The print object should be cleaned with ethanol (96 %) or isopropanol in a cleaning device (e.g. Formlabs Form Wash), or using an ultrasonic bath.

Precautions: When using ultrasound devices, never pour the ethanol (96 %) or isopropanol directly into the ultrasonic bath. Instead, always place it in the recommended container inside the water-filled ultrasonic bath. Use an explosion-proof ultrasonic bath.

Please read carefully the information contained in the instructions for use of the device manufacturer.

1. Clean the print object for 5 min in a reusable cleaning solution using an unheated cleaning device (e.g. Formlabs Form Wash).
2. The print object is then removed from the cleaning bath and sprayed with additional cleaning solution in order to fully rinse off any resin residue.

Tip: The print objects are completely clean when no visible shiny surfaces remain. Resin residue can be removed using a brush soaked in the cleaning solution.

Precautions: The entire cleaning process should not take longer than 5 min, as this could otherwise have a detrimental effect (swelling) on the printed objects.

After cleaning, the print object is dried using compressed air under an extraction unit. If there is liquid resin still adhering to the surface of the object, it can be completely removed by spraying again with cleaning solution and renewed blow drying.

Preparation for post-curing

- Remove support structures. They can be removed using either a cutting disc or side cutters.
- Carbide cutters or diamond grinding stones can be used to make adjustments to the model.

Precautions: A dust mask must also be worn due to potential dust formation while the printed objects are being trimmed/processed.

Post-curing process

The light-curing equipment BEGO Otofash, HiLite Power* and Form Cure* are recommended for post-curing of objects made from VarseoWax Model, as they are optimally suited to the production process.

The final material properties are achieved using a light polymerisation unit with the following performance data: two xenon stroboscopic lamps, flash frequency 10 Hz, light spectrum 300-700 nm (e.g. BEGO Otofash) or one xenon stroboscopic lamp, flash frequency 20 Hz, light spectrum 390-540 nm (e.g. HiLite Power*), or multidirectional LEDs with a flash power 39 W/9.1 W, light spectrum 405 nm (e.g. Form Cure*).

Align/place the models in the equipment with the work surface facing the light source.

An example list of compatible post-curing devices

Light-curing equipment	BEGO Otofash	HiLite Power*	Form Cure*	Note
Flash	2 x 2,000	–	–	Turn object between the exposure cycles
Time	–	2 x 180 sec	2 x 20 min	
Temperature	–	–	60°C	

In order to ensure optimal further processing, please read carefully the information contained in the instructions for use for the light-curing equipment.

The specified times apply only to regularly maintained equipment that can produce the appropriate light intensity.

Tips for further processing of printed and light-cured models made from VarseoWax Model.

1. The models can be insulated to prevent plaster (e.g. Iso 8, BEGO), plastic (e.g. 3D model insulation, Scheu*) and wax (e.g. Isocera, BEGO) from adhering.
2. The models can be duplicated with silicone (e.g. WiroGel®, BEGO) or duplicating gel (e.g. Castogel, BEGO; WiroGel® M, BEGO).
3. It is possible to produce deep drawing copings (Adapta) on stump models.
4. It is possible to produce deep drawing splints on models as a one-off application.

Note: Avoid overheating the foil when deep drawing plastic foils!

5. Repairing printed models made of VarseoWax Model: The models can be repaired with a quick-setting adhesive if a suitable fracture surface for the repair is present.

Note: Refabrication is recommended for large-scale imperfections, cracks or fractures.

In order to ensure optimal processing, please read carefully the information contained in the instructions for use of the different manufacturers.

8. Cleaning in the dental laboratory

Models made from VarseoWax Model can be easily cleaned. Steam cleaning (e.g. with Triton SLA, BEGO) as well as subsequent cleaning with alcohol are possible.

The model can also be cleaned under running water with a toothbrush and soap.

Note: Boiling water (100°C) should not be poured over the model as the model may lose its shape.

9. Disposal

The cured, separated material (base plate, support structure) cannot be re-used. Cured material can be disposed of as domestic waste. Unused resin or ethanol (96 %) used for cleaning or isopropanol with resin residues must be disposed of via the local waste disposal authority or an appropriate hazardous waste collection point, under specification of the safety data sheet.

10. Label symbols

 Manufacturer	 Consult instructions for use
 Date of manufacture	 Use-by date
 LOT Batch code	 Caution
 REF Catalogue number	 Temperature limit
 Keep away from sunlight	

BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG
 Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany
 Tel. +49 421 20 28-0 · Fax +49 421 20 28-100
 E-Mail: info@bego.com · www.bego.com

* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

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