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# BEGO compatibility overview 3D printing system components

## BEGO Varseo XS

### Cleaning

 <p><b>VarseoSmile TriniQ</b> Permanent crowns, bridges and denture teeth</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99%, Ethanol 96% or Inova Printwash)</p>	<p><b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99%)</p>	<p><b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)</p>		
 <p><b>VarseoSmile Crown plus</b> Permanent single crowns, inlays, onlays and veneers</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p><b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99%)</p>	<p><b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/</b> <b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99%)</p>	<p><b>Manuelle cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)</p>
 <p><b>VarseoSmile Temp</b> Temporary crown and bridge restorations, inlays, onlays and veneers</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99%)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99%)</p>	<p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min Trocknung + spraying off (Isopropanol 99%)</p>
 <p><b>VarseoSmile Teeth</b> Denture teeth</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99%)</p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99%) + 3 min reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)</p>	<p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)</p>
 <p><b>VarseoWax CAD/CAST</b> Burnout objects</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)</p>			
 <p><b>VarseoWax Model</b> Dental models</p>	<p><b>Ultrasonic bath:</b> 5 min (Isopropanol 99% or Ethanol 96%)</p>	<p><b>Formlabs Form Wash:</b> 5 min (Isopropanol 99%)</p>	<p><b>Anycubic<sup>1</sup> Wash &amp; Cure Plus<sup>3</sup>:</b> 8 min (Isopropanol 99% or Ethanol 96%)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min reservoir 1 (Isopropanol 99%) + 3 min reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)</p>	

### Light-curing

<p><b>BEGO Otofash:</b> 2 x 2.000 flashes</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters P13 D8 T2.30B0n</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2,5 min exposure parameters P9 D80 T2.00 B0n</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Shining 3D</b> <b>FabCure 2 :</b> 2 x 5 min @ 20 °C</p>		
<p><b>BEGO Otofash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown plus (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox<sup>1</sup> Plus:</b> 2 x 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>
<p><b>Shining 3D</b> <b>FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D</b> <b>FabCure 2 :</b> 2 x 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min Belichtungsparameter: P9 D80 T2.00B0n</p>	<p><b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity Ivl 1 (see preprogrammed resin profile Permanent Crown)</p>	
<p><b>BEGO Otofash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp plus (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox Plus:</b> 2 x 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>
<p><b>Shining 3D</b> <b>FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D</b> <b>FabCure 2 :</b> 2 x 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00B0n</p>	<p><b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity Ivl 1 (see preprogrammed resin profilley CB)</p>	
<p><b>BEGO Otofash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00B0n</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>		
<p><b>Shining 3D</b> <b>FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D</b> <b>FabCure 2 :</b> 2 x 5 min @ 20 °C</p>					
<p><b>BEGO Otofash:</b> 2 x 500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 1 x 2 min exposure parameters P13 D8 T2.00 B0n</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>			
<p><b>Flash or LED light-curing devices, e.g.:</b></p>						
<p><b>BEGO Otofash:</b> 2 x 2.000 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 180 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>Anycubic Wash &amp; Cure Plus<sup>3</sup>:</b> 2 x 20 min</p>	<p><b>Ackuretta CURIE/Varseo Cure<sup>2</sup>:</b> 1 x 5 min exposure parameters: P13 D8 T5.00 B0n</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>3</sup> Compatibility applies to the design status to the serial no. W31126A0405446.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

## BEGO Varseo/Varseo L/Varseo S

### Cleaning

### Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**VarseoSmile Teeth**  
Denture teeth


**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/<sup>2</sup> Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manuelle cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p> <p><b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)</p> <p><b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/<sup>2</sup> Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p> <p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/<sup>2</sup> Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p> <p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/<sup>2</sup> Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p> <p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/<sup>2</sup> Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p> <p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/<sup>2</sup> Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p> <p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/<sup>2</sup> Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p> <p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/<sup>2</sup> Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p> <p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>

<p><b>BEGO Otoflash:</b> 2 × 1.500 flashes</p> <p><b>Shining 3D FabCure:</b> 2 × 15 min</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec</p> <p><b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C</p>	<p><b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C</p> <p><b>Ackuretta CURIE Plus:</b> 2 × 2 min Belichtungsparameter: P9 D80 T2.00B0n</p>	<p><b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C</p> <p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown plus (60 sec – pause – 50 sec) in zone A</p> <p><b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder 12 (25 % Intensität, 405 nm)</p>	<p><b>CUREbox 1 Plus:</b> 2 × 20 min @ 30 °C</p> <p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Ackuretta CURIE/Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30B0n</p> <p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity IV 1 (see preprogrammed resin profile Permanent Crown)</p>
<p><b>BEGO Otoflash:</b> 2 × 1.500 flashes</p> <p><b>Shining 3D FabCure:</b> 2 × 15 min</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec</p> <p><b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C</p>	<p><b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C</p> <p><b>Ackuretta CURIE Plus:</b> 2 × 2 min exposure parameters P9 D80 T2.00B0n</p>	<p><b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C</p> <p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A</p> <p><b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder 12 (25 % Intensität, 405 nm)</p>	<p><b>CUREbox Plus:</b> 2 × 20 min @ 30 °C</p> <p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Ackuretta CURIE/Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30B0n</p> <p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity IV 1 (see preprogrammed resin profiley CB)</p>
<p><b>BEGO Otoflash:</b> 2 × 500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec</p>	<p><b>Ackuretta CURIE/Varseo Cure<sup>2</sup>:</b> 1 × 2 min exposure parameters P13 D8 T2.00 B0n</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>		

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
<sup>2</sup> Both devices are identical in construction.  
<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.  
<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

Ackuretta<sup>1</sup> DENTIQ/Freeshape 120<sup>4</sup>

## Cleaning

## Light-curing

**VarseoSmile TriniQ**

Permanent crowns, bridges and denture teeth

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/</b> <b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manuelle cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
			<p><b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)</p>	<p><b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>

**VarseoSmile Crown Plus**

Permanent single crowns, inlays, onlays and veneers

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p>
			<p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>

**VarseoSmile Temp**

Temporary crown and bridge restorations, inlays, onlays and veneers

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>		
---	--	---	--	--

**VarseoSmile Teeth**

Denture teeth

**VarseoWax CAD/CAST**

Burnout objects

<p><b>Ultrasonic bath:</b> 5 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs Form Wash:</b> 5 min (Isopropanol 99 %)</p>	<p><b>Anycubic<sup>1</sup> Wash &amp; Cure Plus<sup>3</sup>:</b> 8 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
---	--	---	---

**VarseoWax Model**

Dental models

<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown Plus (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox<sup>1</sup> Plus:</b> 2 x 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>
<p><b>Shining 3D FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min Belichtungparameter: P9 D80 T2.00B0n</p>	<p><b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>		<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity lvl 1 (see preprogrammed resin profile Permanent Crown)</p>

<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox Plus:</b> 2 x 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>
<p><b>Shining 3D FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00B0n</p>	<p><b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>		<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity lvl 1 (see preprogrammed resin profiley CB)</p>

<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00B0n</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>		
<p><b>Shining 3D FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C</p>					

<b>Flash or LED light-curing devices, e.g.:</b>						
<p><b>BEGO Otoflash:</b> 2 x 2.000 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 180 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>Anycubic Wash &amp; Cure Plus<sup>3</sup>:</b> 2 x 20 min</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 1 x 5 min exposure parameters: P13 D8 T5.00 B0n</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
<sup>2</sup> Both devices are identical in construction.  
<sup>3</sup> Compatibility applies to the design status to the serial no. W31126A0405446.  
<sup>4</sup> Ackuretta Dentiq & Freeshape 120 can only be used with anodized aluminium build platform (for VarseoSmile Crown Plus and VarseoSmile Temp).  
<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.  
<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

Ackuretta<sup>1</sup> SOL<sup>5</sup>

## Cleaning

## Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 %, Ethanol 96 % or Inova Printwash)	<b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
--	--	---

<b>BEGO Otoflash:</b> 2 x 2.000 flashes	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters P13 D8 T2.30BOn	<b>Ackuretta CURIE Plus:</b> 2 x 2,5 min exposure parameters P9 D80 T2.00 BOn	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C
--	---	--	--	--

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/ Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manuelle cleaning with tensides:</b> Details see last page	<b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
		<b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)		<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

<b>BEGO Otoflash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown plus (60 sec – pause – 50 sec) in zone A	<b>CUREbox<sup>1</sup> Plus:</b> 2 x 20 min @ 30 °C	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min Belichtungsparameter: P9 D80 T2.00BOn		<b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder <sup>12</sup> (25% Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity lvl 1 (see preprogrammed resin profile Permanent Crown)

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix Veriwash + Veriwhirl/ Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manual cleaning with tensides:</b> Details see last page	<b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)
				<b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
			<b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)	

<b>BEGO Otoflash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer):</b> 2 x 90 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A	<b>CUREbox Plus:</b> 2 x 20 min @ 30 °C	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00BOn		<b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder <sup>12</sup> (25% Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity lvl 1 (see preprogrammed resin profile CB)

**VarseoSmile Teeth**  
Denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
---	---	---

<b>BEGO Otoflash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer):</b> 2 x 90 sec	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00BOn	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C			

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
<sup>2</sup> Both devices are identical in construction.  
<sup>5</sup> Ackuretta SOL can only be used with anodized aluminium build platform (Small, Medium and Large).  
<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.  
<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

Anycubic<sup>1</sup> Photon Mono X<sup>6</sup>

## Cleaning

## Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<b>Ultrasonic bath:</b> 5 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 5 min (Isopropanol 99 %)	<b>Anycubic<sup>1</sup> Wash &amp; Cure Plus<sup>3</sup>:</b> 8 min (Isopropanol 99 % or Ethanol 96 %)	<b>Rapid Shape RS Wash:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
--	---	--	--

<b>Flash or LED light-curing devices, e.g.:</b>					
<b>BEGO Otoflash:</b> 2 × 2.000 flashes	<b>HiLite Power (Kulzer):</b> 2 × 180 sec	<b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C	<b>Anycubic Wash &amp; Cure Plus<sup>3</sup>:</b> 2 × 20 min	<b>Ackuretta CURIE/Varseo Cure<sup>2</sup>:</b> 1 × 5 min exposure parameters: P13 D8 T5.00 B0n	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
<sup>3</sup> Compatibility applies to the design status to the serial no. W31126A0405446.  
<sup>6</sup> Compatibility applies to the design status to the serial no. P02123C0306508.  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

## Asiga<sup>1</sup> Max UV/Max 405

### Cleaning

 <p><b>VarseoSmile TriniQ</b> Permanent crowns, bridges and denture teeth</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 %, Ethanol 96 % or Inova Printwash)</p>	<p><b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>		
 <p><b>VarseoSmile Crown plus</b> Permanent single crowns, inlays, onlays and veneers</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/</b> <b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manuelle cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
 <p><b>VarseoSmile Temp</b> Temporary crown and bridge restorations, inlays, onlays and veneers</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p>
 <p><b>VarseoSmile Teeth</b> Denture teeth</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
 <p><b>VarseoWax CAD/CAST</b> Burnout objects</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>			
 <p><b>VarseoWax Model</b> Dental models</p>	<p><b>Ultrasonic bath:</b> 5 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs Form Wash:</b> 5 min (Isopropanol 99 %)</p>	<p><b>Anycubic<sup>1</sup> Wash &amp; Cure Plus<sup>3</sup>:</b> 8 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	

### Light-curing

<p><b>BEGO Otoflash:</b> 2 x 2.000 flashes</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters P13 D8 T2.30B0n</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2,5 min exposure parameters P9 D80 T2.00 B0n</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Shining 3D</b> <b>FabCure 2 :</b> 2 x 5 min @ 20 °C</p>		
<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown plus (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox 1 Plus:</b> 2 x 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>
<p><b>Shining 3D</b> <b>FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D</b> <b>FabCure 2 :</b> 2 x 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min Belichtungsparameter: P9 D80 T2.00B0n</p>	<p><b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder 12 (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity Ivl 1 (see preprogrammed resin profile Permanent Crown)</p>	
<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox Plus:</b> 2 x 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>
<p><b>Shining 3D</b> <b>FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D</b> <b>FabCure 2 :</b> 2 x 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00B0n</p>	<p><b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder 12 (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity Ivl 1 (see preprogrammed resin profiley CB)</p>	
<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00B0n</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30B0n</p>		
<p><b>Shining 3D</b> <b>FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D</b> <b>FabCure 2 :</b> 2 x 5 min @ 20 °C</p>					
<p><b>BEGO Otoflash:</b> 2 x 500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 1 x 2 min exposure parameters P13 D8 T2.00 B0n</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>			
Flash or LED light-curing devices, e.g.:						
<p><b>BEGO Otoflash:</b> 2 x 2.000 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 180 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>Anycubic Wash &amp; Cure Plus<sup>3</sup>:</b> 2 x 20 min</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 1 x 5 min exposure parameters: P13 D8 T5.00 B0n</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>3</sup> Compatibility applies to the design status to the serial no. W31126A0405446.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

Asiga<sup>1</sup> Pro 4K80<sup>7</sup>

## Cleaning

## Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<b>Ultrasonic bath:</b> 5 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 5 min (Isopropanol 99 %)	<b>Anycubic<sup>1</sup> Wash &amp; Cure Plus<sup>3</sup>:</b> 8 min (Isopropanol 99 % or Ethanol 96 %)	<b>Rapid Shape RS Wash:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
--	---	--	---

<b>BEGO Otoflash:</b> 2 × 2.000 flashes	<b>HiLite Power (Kulzer):</b> 2 × 180 sec	<b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C	<b>Anycubic Wash &amp; Cure Plus<sup>3</sup>:</b> 2 × 20 min	<b>Ackuretta CURIE/Varseo Cure<sup>2</sup>:</b> 1 × 5 min exposure parameters: P13 D8 T5.00 B0n	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)
--	--	--	---	--	--

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
<sup>3</sup> Compatibility applies to the design status to the serial no. W31126A0405446.  
<sup>7</sup> Printer must be operated in 4K mode.  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

DEKEMA<sup>1</sup> trix print<sup>2</sup>

## Cleaning

## Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/ Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manuelle cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
			<p><b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)</p>	<p><b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p>
			<p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>		
---	---	---	--	--

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown plus (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox<sup>1</sup> Plus:</b> 2 x 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
<p><b>Shining 3D FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D FabCure 2:</b> 2 x 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min Belichtungsparameter: P9 D80 T2.00BOn</p>		<p><b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity Ivl 1 (see preprogrammed resin profile Permanent Crown)</p>

<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox Plus:</b> 2 x 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
<p><b>Shining 3D FabCure:</b> 2 x 15 min</p>	<p><b>Shining 3D FabCure 2:</b> 2 x 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00BOn</p>		<p><b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity Ivl 1 (see preprogrammed resin profiley CB)</p>

<p><b>BEGO Otoflash:</b> 2 x 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p>	<p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00BOn</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
--	---	---	---	--

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
<sup>2</sup> Both devices are identical in construction.  
<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.  
<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

Formlabs<sup>1</sup> Form 2<sup>9</sup> (stainless steel build platform necessary) 

Formlabs Form 3B/3B+ (stainless steel build platform necessary) 

## Cleaning

## Light-curing

## Cleaning

## Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

**Formlabs Form Wash:**  
3 min (Isopropanol 99%)

**Formlabs Form Cure:**  
2 × 20 min @ 60 °C  
  
**Formlabs Fast Cure<sup>10</sup>:**  
2 × 2.5 min with intensity lvl 1  
(see preprogrammed resin profile Permanent Crown)

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**Formlabs Form Wash:**  
3 min (Isopropanol 99%)

**Formlabs Form Cure:**  
2 × 20 min @ 60 °C  
  
**Formlabs Fast Cure<sup>10</sup>:**  
2 × 2.5 min with intensity lvl 1  
(see preprogrammed resin profile Temporary CB)

**Formlabs Form Wash:**  
3 min (Isopropanol 99%)

**Formlabs Form Cure:**  
2 × 20 min @ 60 °C  
  
**Formlabs Fast Cure<sup>10</sup>:**  
2 × 2.5 min with intensity lvl 1  
(see preprogrammed resin profile Temporary CB)

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>9</sup> VarseoSmile Crown plus is distributed by Formlabs as Permanent Crown + VarseoSmile Temp as Temporary CB.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>13</sup> Bleach shade only available for Form 3B/+

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

## Ivoclar PrograPrint PR5

### Cleaning

### Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/ Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manuelle cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
			<p><b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)</p>	<p><b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p>
			<p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>

<p><b>BEGO Otoflash:</b> 2 × 1.500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown plus (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox<sup>1</sup> Plus:</b> 2 × 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
<p><b>Shining 3D FabCure:</b> 2 × 15 min</p>	<p><b>Shining 3D FabCure 2:</b> 2 × 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 × 2 min Belichtungsparameter: P9 D80 T2.00BOn</p>	<p><b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity Ivl 1 (see preprogrammed resin profile Permanent Crown)</p>	
<p><b>BEGO Otoflash:</b> 2 × 1.500 flashes</p>	<p><b>HiLite Power (Kulzer):</b> 2 × 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox Plus:</b> 2 × 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
<p><b>Shining 3D FabCure:</b> 2 × 15 min</p>	<p><b>Shining 3D FabCure 2:</b> 2 × 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 × 2 min exposure parameters P9 D80 T2.00BOn</p>	<p><b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity Ivl 1 (see preprogrammed resin profilel CB)</p>	

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

## Microlay<sup>1</sup> Versus 385

### Cleaning

### Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/</b> <b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manuelle cleaning with tensides:</b> <i>Details see last page</i>	<b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
			<b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)	<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

<b>BEGO Otofash:</b> 2 × 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec	<b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown <sup>plus</sup> (60 sec – pause – 50 sec) in zone A	<b>CUREbox<sup>1</sup> Plus:</b> 2 × 20 min @ 30 °C	<b>Ackuretta CURIE / Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30B0n
<b>Shining 3D FabCure:</b> 2 × 15 min	<b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 × 2 min Belichtungsparameter: P9 D80 T2.00B0n		<b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder <sup>12</sup> (25 % Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity Ivl 1 (see preprogrammed resin profile Permanent Crown)

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manual cleaning with tensides:</b> <i>Details see last page</i>	<b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)
			<b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)	<b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

<b>BEGO Otofash:</b> 2 × 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec	<b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A	<b>CUREbox Plus:</b> 2 × 20 min @ 30 °C	<b>Ackuretta CURIE / Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30B0n
<b>Shining 3D FabCure:</b> 2 × 15 min	<b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 × 2 min exposure parameters P9 D80 T2.00B0n		<b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder <sup>12</sup> (25 % Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity Ivl 1 (see preprogrammed resin profile CB)

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

<b>Ultrasonic bath:</b> 5 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 5 min (Isopropanol 99 %)	<b>Anycubic<sup>1</sup> Wash &amp; Cure Plus<sup>3</sup>:</b> 8 min (Isopropanol 99 % or Ethanol 96 %)	<b>Rapid Shape RS Wash:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
--	---	--	---

<b>Flash or LED light-curing devices, e.g.:</b>						
<b>BEGO Otofash:</b> 2 × 2.000 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 180 sec	<b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C	<b>Anycubic Wash &amp; Cure Plus<sup>3</sup>:</b> 2 × 20 min	<b>Ackuretta CURIE / Varseo Cure<sup>2</sup>:</b> 1 × 5 min exposure parameters: P13 D8 T5.00 B0n	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	

**VarseoWax Model**  
Dental models

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>3</sup> Compatibility applies to the design status to the serial no. W31126A0405446.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

Planmeca Creo C5 <sup>13</sup>

## Cleaning

## Light-curing

 <p><b>VarseoSmile TriniQ</b> Permanent crowns, bridges and denture teeth</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs <sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix <sup>1</sup> Veriwash + Veriwhirl/ Ackuretta <sup>1</sup> Cleani <sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manuelle cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay <sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>Ivoclar <sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)</p> <p><b>Rapid Shape <sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>BEGO Otofash:</b> 2 x 1.500 flashes</p> <p><b>HiLite Power (Kulzer <sup>1</sup>):</b> 2 x 90 sec</p> <p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p> <p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p> <p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown <sup>plus</sup> (60 sec – pause – 50 sec) in zone A</p> <p><b>CUREbox <sup>1</sup> Plus:</b> 2 x 20 min @ 30 °C</p> <p><b>Ackuretta CURIE/ Varseo Cure <sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
 <p><b>VarseoSmile Crown <sup>plus</sup></b> Permanent single crowns, inlays, onlays and veneers</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani <sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p>	<p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>BEGO Otofash:</b> 2 x 1.500 flashes</p> <p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p> <p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p> <p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p> <p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp <sup>plus</sup> (60 sec – pause – 50 sec) in zone A</p> <p><b>CUREbox Plus:</b> 2 x 20 min @ 30 °C</p> <p><b>Ackuretta CURIE/ Varseo Cure <sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
 <p><b>VarseoSmile Temp</b> Temporary crown and bridge restorations, inlays, onlays and veneers</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani <sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p>	<p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>BEGO Otofash:</b> 2 x 1.500 flashes</p> <p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p> <p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p> <p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p> <p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp <sup>plus</sup> (60 sec – pause – 50 sec) in zone A</p> <p><b>CUREbox Plus:</b> 2 x 20 min @ 30 °C</p> <p><b>Ackuretta CURIE/ Varseo Cure <sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
 <p><b>VarseoSmile Teeth</b> Denture teeth</p>	<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> <p><b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani <sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p> <p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p> <p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>	<p><b>BEGO Otofash:</b> 2 x 1.500 flashes</p> <p><b>HiLite Power (Kulzer):</b> 2 x 90 sec</p> <p><b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C</p> <p><b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C</p> <p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp <sup>plus</sup> (60 sec – pause – 50 sec) in zone A</p> <p><b>CUREbox Plus:</b> 2 x 20 min @ 30 °C</p> <p><b>Ackuretta CURIE/ Varseo Cure <sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
 <p><b>VarseoWax CAD/CAST</b> Burnout objects</p>			
 <p><b>VarseoWax Model</b> Dental models</p>			

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

<sup>13</sup> Creo C5 can only be used with Glass Build Platform.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

## Original Prusa<sup>1</sup> Medical One

### Cleaning

### Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/</b> <b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manuelle cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>
			<p><b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)</p>	<p><b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>

<p><b>BEGO Otoflash:</b> 2 × 1.500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown<sup>plus</sup> (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox<sup>1</sup> Plus:</b> 2 × 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
<p><b>Shining 3D FabCure:</b> 2 × 15 min</p>	<p><b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 × 2 min Belichtungsparameter: P9 D80 T2.00BOn</p>		<p><b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity Ivl 1 (see preprogrammed resin profile Permanent Crown)</p>

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>Manual cleaning with tensides:</b> <i>Details see last page</i></p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)</p>
			<p><b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)</p>	<p><b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>

<p><b>BEGO Otoflash:</b> 2 × 1.500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec</p>	<p><b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C</p>	<p><b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C</p>	<p><b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A</p>	<p><b>CUREbox Plus:</b> 2 × 20 min @ 30 °C</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
<p><b>Shining 3D FabCure:</b> 2 × 15 min</p>	<p><b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 × 2 min exposure parameters P9 D80 T2.00BOn</p>		<p><b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder<sup>12</sup> (25 % Intensität, 405 nm)</p>	<p><b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)</p>	<p><b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity Ivl 1 (see preprogrammed resin profiley CB)</p>

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

<p><b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p><b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)</p>	<p><b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p>		
---	--	---	--	--

<p><b>BEGO Otoflash:</b> 2 × 1.500 flashes</p>	<p><b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec</p>	<p><b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C</p>	<p><b>Ackuretta CURIE Plus:</b> 2 × 2 min exposure parameters P9 D80 T2.00BOn</p>	<p><b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30BOn</p>
<p><b>Shining 3D FabCure:</b> 2 × 15 min</p>	<p><b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C</p>			

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

## Rapid Shape D10+<sup>15</sup>/D20 II/D20+/D30 II/D30+/D50+

### Cleaning

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 %, Ethanol 96 % or Inova Printwash)	<b>Ackuretta<sup>1</sup> Cleanj<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
---	---	--

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl:</b> <b>Ackuretta<sup>1</sup> Cleanj<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manuelle cleaning with tensides:</b> Details see last page	<b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
			<b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)	<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

**VarseoSmile Crown Plus**  
Permanent single crowns, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix Veriwash + Veriwhirl:</b> <b>Ackuretta<sup>1</sup> Cleanj<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manual cleaning with tensides:</b> Details see last page	<b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)
			<b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)	<b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**VarseoSmile Teeth**  
Denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
--	--

**VarseoWax CAD/CAST**  
Burnout objects

<b>Ultrasonic bath:</b> 5 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 5 min (Isopropanol 99 %)	<b>Anycubic<sup>1</sup> Wash &amp; Cure Plus<sup>3</sup>:</b> 8 min (Isopropanol 99 % or Ethanol 96 %)	<b>Rapid Shape RS Wash:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
--	---	--	--

**VarseoWax Model**  
Dental models

### Light-curing

<b>BEGO Otofash:</b> 2 x 2.000 flashes	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2.5 min exposure parameters P13 D8 T2.30B0n	<b>Ackuretta CURIE Plus:</b> 2 x 2.5 min exposure parameters P9 D80 T2.00 B0n	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C
---	---	--	--	---

<b>BEGO Otofash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown <sup>Plus</sup> (60 sec – pause – 50 sec) in zone A	<b>CUREbox<sup>1</sup> Plus:</b> 2 x 20 min @ 30 °C	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2.5 min exposure parameters: P13 D8 T2.30B0n
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min Belichtungsparameter: P9 D80 T2.00B0n	<b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder <sup>12</sup> (25 % Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2.5 min with intensity lvl 1 (see preprogrammed resin profile Permanent Crown)	

<b>BEGO Otofash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer):</b> 2 x 90 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A	<b>CUREbox Plus:</b> 2 x 20 min @ 30 °C	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2.5 min exposure parameters: P13 D8 T2.30B0n
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2x 2 min exposure parameters P9 D80 T2.00B0n	<b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder <sup>12</sup> (25 % Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2.5 min with intensity lvl 1 (see preprogrammed resin profilel CB)	

<b>BEGO Otofash:</b> 2 x 500 flashes	<b>HiLite Power (Kulzer):</b> 2 x 90 sec	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 1 x 2 min exposure parameters P13 D8 T2.00 B0n	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)
---	---	--	--

<b>BEGO Otofash:</b> 2 x 2.000 flashes	<b>HiLite Power (Kulzer):</b> 2 x 180 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>Anycubic Wash &amp; Cure Plus<sup>3</sup>:</b> 2 x 20 min	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 1 x 5 min exposure parameters: P13 D8 T5.00 B0n	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)
---	--	---	---	---	--

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>3</sup> Compatibility applies to the design status to the serial no. W31126A0405446.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

<sup>15</sup> VarseoSmile TriniQ is not compatible with RapidShape D10+  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.  
 Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.  
 Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

Shining 3D AccuFab D1S/AccuFab L4D <sup>14 16</sup>

## Cleaning

## Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 %, Ethanol 96 % or Inova Printwash)	<b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
--	--	---

<b>BEGO Otofash:</b> 2 × 2.000 flashes	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters P13 D8 T2.30BOn	<b>Ackuretta CURIE Plus:</b> 2 × 2,5 min exposure parameters P9 D80 T2.00 BOn	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C
---	---	--	--	--

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/ Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manuelle cleaning with tensides:</b> <i>Details see last page</i>	<b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
			<b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)	<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

<b>BEGO Otofash:</b> 2 × 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 × 90 sec	<b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown plus (60 sec – pause – 50 sec) in zone A	<b>CUREbox<sup>1</sup> Plus:</b> 2 × 20 min @ 30 °C	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 × 15 min	<b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 × 2 min Belichtungsparameter: P9 D80 T2.00BOn	<b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder <sup>12</sup> (25 % Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity Ivl 1 (see preprogrammed resin profile Permanent Crown)	

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manual cleaning with tensides:</b> <i>Details see last page</i>	<b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)
			<b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)	<b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

<b>BEGO Otofash:</b> 2 × 1.500 flashes	<b>HiLite Power (Kulzer):</b> 2 × 90 sec	<b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A	<b>CUREbox Plus:</b> 2 × 20 min @ 30 °C	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 × 15 min	<b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 × 2 min exposure parameters P9 D80 T2.00BOn	<b>Ivoclar PrograPrint Cure:</b> 2 × 2 min on the PrograPrint Object Holder <sup>12</sup> (25 % Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 × 2,5 min with intensity Ivl 1 (see preprogrammed resin profiley CB)	

**VarseoSmile Teeth**  
Denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
---	---	---

<b>BEGO Otofash:</b> 2 × 1.500 flashes	<b>HiLite Power (Kulzer):</b> 2 × 90 sec	<b>SprintRay ProCure:</b> 2 × 20 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 × 2 min exposure parameters P9 D80 T2.00BOn	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 × 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 × 15 min	<b>Shining 3D FabCure 2 :</b> 2 × 5 min @ 20 °C			

**VarseoWax CAD/CAST**  
Burnout objects

--	--	--

--	--	--	--	--

**VarseoWax Model**  
Dental models

--	--	--

--	--	--	--	--

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

<sup>14</sup> AccuFab D1S and AccuFab L4D can only be used with Ceramic Build Platform.

<sup>16</sup> VarseoSmile TriniQ is only validated for use with AccuFab L4D

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

SprintRay<sup>1</sup> Pro 95/Pro 95 S/Pro 55/Pro 55 S

## Cleaning

## Light-curing

**VarseoSmile TriniQ**

Permanent crowns, bridges and denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99% or Ethanol 96%)	<b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99%)	<b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl/</b> <b>Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99%)	<b>Manuelle cleaning with tensides:</b> <i>Details see last page</i>	<b>SprintRay<sup>1</sup> ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)
			<b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + brush cleaning + spray off (Isopropanol 99%)	<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)

**VarseoSmile Crown plus**

Permanent single crowns, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99% or Ethanol 96%)	<b>Formlabs Form Wash:</b> 3 min (Isopropanol 99%)	<b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99%)	<b>Manual cleaning with tensides:</b> <i>Details see last page</i>	<b>SprintRay ProWash S/ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min Trocknung + spraying off (Isopropanol 99%)
			<b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + brush cleaning + spraying off (Isopropanol 99%)	<b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99%) + 3 min Reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)

**VarseoSmile Temp**

Temporary crown and bridge restorations, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99% or Ethanol 96%)	<b>Whip Mix Veriwash + Veriwhirl/</b> <b>Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99%)	<b>SprintRay ProWash S/ProWash/Dry:</b> 4 min reservoir 1 (Isopropanol 99%) + 3 min reservoir 2 (Isopropanol 99%) + 3 min drying + spraying off (Isopropanol 99%)		
---	--	--	--	--

**VarseoSmile Teeth**

Denture teeth

**VarseoWax CAD/CAST**

Burnout objects

**VarseoWax Model**

Dental models

<b>BEGO Otoflash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown plus (60 sec – pause – 50 sec) in zone A	<b>CUREbox<sup>1</sup> Plus:</b> 2 x 20 min @ 30 °C	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2:</b> 2 x 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min Belichtungsparameter: P9 D80 T2.00BOn	<b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder <sup>12</sup> (25% Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100% Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity I1 (see preprogrammed resin profile Permanent Crown)	

<b>BEGO Otoflash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A	<b>CUREbox Plus:</b> 2 x 20 min @ 30 °C	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2:</b> 2 x 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00BOn	<b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder <sup>12</sup> (25% Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100% Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity I1 (see preprogrammed resin profile CB)	

<b>BEGO Otoflash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00BOn	<b>Ackuretta CURIE/ Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn		
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2:</b> 2 x 5 min @ 20 °C					

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
<sup>2</sup> Both devices are identical in construction.  
<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.  
<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

Whip Mix<sup>1</sup> VeriBuild/VeriEko<sup>11</sup>

## Cleaning

## Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs<sup>1</sup> Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix<sup>1</sup> Veriwash + Veriwhirl / Ackuretta<sup>1</sup> Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manuelle cleaning with tensides:</b> <i>Details see last page</i>	<b>SprintRay<sup>1</sup> ProWash S / ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
			<b>Ivoclar<sup>1</sup> PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spray off (Isopropanol 99 %)	<b>Rapid Shape<sup>1</sup> RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

<b>BEGO Otofflash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer<sup>1</sup>):</b> 2 x 90 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Crown <sup>plus</sup> (60 sec – pause – 50 sec) in zone A	<b>CUREbox<sup>1</sup> Plus:</b> 2 x 20 min @ 30 °C	<b>Ackuretta CURIE / Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min Belichtungsparameter: P9 D80 T2.00BOn	<b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder <sup>12</sup> (25 % Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity Ivl 1 (see preprogrammed resin profile Permanent Crown)	

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

<b>Ultrasonic bath:</b> 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 3 min (Isopropanol 99 %)	<b>Whip Mix Veriwash + Veriwhirl / Ackuretta Cleani<sup>2</sup>:</b> 3 min + 3 min (Isopropanol 99 %)	<b>Manual cleaning with tensides:</b> <i>Details see last page</i>	<b>SprintRay ProWash S / ProWash/Dry:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min Trocknung + spraying off (Isopropanol 99 %)
			<b>Ivoclar PrograPrint Clean:</b> 3 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)	<b>Rapid Shape RS Wash:</b> 4 min Reservoir 1 (Isopropanol 99 %) + 3 min Reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

<b>BEGO Otofflash:</b> 2 x 1.500 flashes	<b>HiLite Power (Kulzer):</b> 2 x 90 sec	<b>Formlabs Form Cure:</b> 2 x 20 min @ 60 °C	<b>SprintRay ProCure:</b> 2 x 20 min @ 20 °C	<b>SprintRay ProCure 2:</b> Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec) in zone A	<b>CUREbox Plus:</b> 2 x 20 min @ 30 °C	<b>Ackuretta CURIE / Varseo Cure<sup>2</sup>:</b> 2 x 2,5 min exposure parameters: P13 D8 T2.30BOn
<b>Shining 3D FabCure:</b> 2 x 15 min	<b>Shining 3D FabCure 2 :</b> 2 x 5 min @ 20 °C	<b>Ackuretta CURIE Plus:</b> 2 x 2 min exposure parameters P9 D80 T2.00BOn	<b>Ivoclar PrograPrint Cure:</b> 2 x 2 min on the PrograPrint Object Holder <sup>12</sup> (25 % Intensität, 405 nm)	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)	<b>Formlabs Fast Cure<sup>10</sup>:</b> 2 x 2,5 min with intensity Ivl 1 (see preprogrammed resin profily CB)	

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

<sup>2</sup> Both devices are identical in construction.

<sup>10</sup> Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

<sup>11</sup> Whip Mix VeriBuild & VeriEKO can only be used with anodized aluminium build platform (VeriEKO anodized aluminium build platform in Small, Medium and Large).

<sup>12</sup> PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# BEGO compatibility overview 3D printing system components

## DLP printers in general (with 385 – 405 nm wave length)

### Cleaning

### Light-curing

**VarseoSmile TriniQ**  
Permanent crowns, bridges and denture teeth

**VarseoSmile Crown plus**  
Permanent single crowns, inlays, onlays and veneers

**VarseoSmile Temp**  
Temporary crown and bridge restorations, inlays, onlays and veneers

**VarseoSmile Teeth**  
Denture teeth

**VarseoWax CAD/CAST**  
Burnout objects

**VarseoWax Model**  
Dental models

<b>Ultrasonic bath:</b> 5 min (Isopropanol 99 % or Ethanol 96 %)	<b>Formlabs Form Wash:</b> 5 min (Isopropanol 99 %)	<b>Anycubic<sup>1</sup> Wash &amp; Cure Plus<sup>3</sup>:</b> 8 min (Isopropanol 99 % or Ethanol 96 %)	<b>Rapid Shape RS Wash:</b> 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)
--	---	--	--

<b>BEGO Otoflash:</b> 2 × 2.000 flashes	<b>HiLite Power (Kulzer):</b> 2 × 180 sec	<b>Formlabs Form Cure:</b> 2 × 20 min @ 60 °C	<b>Anycubic Wash &amp; Cure Plus<sup>3</sup>:</b> 2 × 20 min	<b>Ackuretta CURIE/Varseo Cure<sup>2</sup>:</b> 1 × 5 min exposure parameters: P13 D8 T5.00 B0n	<b>Rapid Shape RS Cure:</b> 15 min using lower and upper wavelength at 100 % Power (without Vakuum)
--	--	--	---	--	--

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
<sup>2</sup> Compatibility applies to the design status to the serial no. W31126A0405446.  
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

# Manual cleaning with Tensides

This cleaning method is valid for VarseoSmile Crown<sup>plus</sup> and VarseoSmile Temp


## Necessary tools, equipment and materials

- InovaPrint wash (hp-dent<sup>1</sup>) general purpose 3D print cleaner
- Tap water
- 1-propanol (70 Vol.-%)
- Toothbrush
- Brush
- Instrument for holding the printed objects (e.g. artery clamp)
- Absorbent pad (e.g. paper towels) or tub
- 2 cups for 1-propanol (70 Vol.-%)
- 1 cup for cleaning solution
- Compressed air with trigger/splash guard
- Personal protective equipment: protective gloves and goggles

## Cleaning process

Carefully remove the printed object from the platform with the help of a spatula.

Clean the printed object by using a toothbrush and/or brush in two steps:

### Step 1:

**Pre-wash with reusable cleaning solution: 5 % InovaPrint wash (hp-dent) + 95 % tap water.**

- 1.1 Vigorously swirl object in cleaning solution for 15 s using pliers or an arterial clamp. Ensure not to damage the printed object.

- 1.2 Remove excess resin by using a toothbrush and cleaning solution. In addition, a brush can be used to clean the inside of crowns. Shortly swirl object in cleaning solution as needed.

**Note:** Clean toothbrush and brush regularly with 1-Propanol when covered with resin. Dry before reusing toothbrush or brush again for cleaning the printed object.

- 1.3 Use compressed air under a fume hood with splash guard to remove the cleaning solution from the printed object's surface.

- 1.4 Repeat step 1.2 and 1.3 until just a thin layer of resin remains on the surface of the print.

### Step 2:

**Finish with fresh 1-Propanol (70 Vol.-%).**

- 2.1 Vigorously swirl printed object in 1-Propanol for 5 s and dry object immediately with compressed air.


- 2.2 Check for remaining resin (shiny spots). If surface of printed object is matte, remove support structure as described in the following step (step 2.3). Shiny spots can be removed by quickly brushing the surface with a 1-Propanol soaked brush. Immediately dry object with compressed air.

- 2.3 Remove the support structure with the help of a cutting wheel or side cutters. Ensure not to deform the printed object.

- 2.4 Quickly clean occlusal surface with a 1-Propanol soaked brush to remove any excess resin. Immediately dry printed object with compressed air.

**Note:** Contact time of printed object with 1-Propanol needs to be reduced to a minimum to avoid the formation of white spots on the surface.

- 2.5 Post-cure the dental objects by using a validated post-curing unit.

<sup>1</sup> This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.